

inspection schedule

Heidelberg Upright Stander Jumbo(hip-knee-ankle-foot orthosis)

operator

Product

serial number

Scheduled maintenance intervals

--

--

12 months

Item	Assembly		Settings & functions		Damage & deformation		Screw connections	
			without defects	with defects	without defects	with defects	without defects	with defects
1	Frame	Base frame						
		Transport castors and parking brake (optional)						
		Feet						
		Wall and furniture protection						
		Clamping and fastening elements						
		Multi base frame (optional)						
		Centre column (table column)						
2	Foot area	Footplate						
		Foot shells (optional)*						
		Foot straps (optional)*						
		Clamping and fastening elements						
3	Knee area*	Knee pads Standard/Anatomical						
		Knee pad holders						
		Knee pad padding and covers						
		Clamping and fastening elements						
4	Pelvic area	Pelvic pads*						
		- Brackets						
		- Upholstery and covers						
		- Fastening/clamping elements						
		Swivel arm sink frame						
		- Brackets						
		- Upholstery and covers						
		- Fastening/clamping elements						
		Pillar						
		- Brackets / mounts						
		- Fastening/clamping elements						

(*) – not for the Heidelberg Jumbo for hip-knee-ankle-foot orthosis

Heidelberg Upright Stander Jumbo(hip-knee-ankle-foot orthosis)

Item	Assembly		Settings & functions		Damage & deformation		Screw connections	
			without defects	with defects	without defects	with defects	without defects	with defects
5	Thoracic area	Chest pad (sternum; curved shape, made to measure)						
		- Brackets						
		- Pads and covers						
		- Fastening/clamping elements						
6	Back area	Back pad						
		- Brackets						
		- Pads and covers						
		- Fastening/clamping elements						
		Headrests (optional)						
		- Brackets						
		- Pads and covers						
		- Fastening/clamping elements						
7	Table area	Table area / armrests						
		- Brackets						
		- Fastening eyelets (optional)**						
		- Table mount with depth adjustment						
		- Wooden parts						
		- Fastening/clamping elements						

(**) – only for the Heidelberg Jumbo standing frame for standing orthoses

Notes for any repairs and further maintenance:

Inspection carried out on

Inspection carried out by

Signature