

Operator

Product

--

Serial number

--

Maintenance interval

12 months

Pos.	Assembly		Settings & functions		Damage & deformation		Bolted joints	
			no defects	defects found	no defects	defects found	no defects	defects found
1	Frame	Standard frame, straight						
		Frame, abducted						
		Long frame						
		Tip assist						
		Parking brakes/knee lever brakes						
		Castor wheels with wheel forks and castor stems						
		Direction lock						
		Drive wheels with drum brake carrier and hand brake levers						
		Quick-release axle						
		Spokeguard						
		Wheel guard						
		Seat tilting via gas spring						
		Wheelchair-tiedown and occupant-restraint system						
		scala-/viamobil attachment						
		scalamobil anti-tip						
2	Seat	Seat plate						
		Trapezoid adapter						
		Mounting for footplate						
3	Back	Back adjustment, mechanical						
		Back adjustment via gas spring						
		Push bar, adjustable for height and angle						
		Push handles						
		Back mounting system with a split, fixed angle						
		Back mounting system with a split, adjustable angle						
		Standard armrests						
		Comfort armrests						

Pos.	Assembly		Settings & functions		Damage & deformation		Bolted joints	
			no defects	defects found	no defects	defects found	no defects	defects found
4	Feet	Footplate						
		Split footplate						
		Adjustable foot angle						
		Folding mechanism						
		Heel loop						
		Calf supports						
		Split calf supports						

Notes on repairs and maintenance

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There is no handwriting or other markings on the paper.

Inspection completed on

[illegible]

Inspection completed by

Signature

Pos.	Assembly		Settings & functions		Damage & deformation		Bolted joints	
			no defects	defects found	no defects	defects found	no defects	defects found

Notes on repairs and maintenance

Inspection completed on

Inspection completed by

Signature