

# Inspection Schedule

## Heidelberg Upright Stander Benni

operator


Product

Serial number

Scheduled maintenance intervals

--

--

<b>12 months</b>
------------------

Item	Assembly		Settings & functions		Damage & deformation		Screw connections	
			without defects	with defects	without defects	with defects	without defects	with defects
1	Frame	base frame						
		Transport castors and parking brake (optional)						
		Feet						
		Wall and furniture protection (optional)						
		Rubber damper						
		Blocker and carrier plate						
		Clamping and fastening elements						
		Dynamic centre column						
2	Foot area	Footboard						
		Footrests						
		Foot strap (optional)						
		Clamping and fastening elements						
3	Knee area	Knee braces Standard/Anatomical						
		Brackets for knee pads						
		Padding and covers for knee pads						
		Clamping and fastening elements						
4	Pelvic area	Pelvic pad						
		- Brackets						
		- Upholstery and covers						
		- Fastening/clamping elements						
		Swivel arm frame						
		- Frame						
		- Brackets						
		- Fastening/clamping elements						

## Heidelberg Upright Stander Benni

Item	Assembly		Settings & functions		Damage & deformation		Screw connections	
			without defects	with defects	without defects	with defects	without defects	with defects
		Dynamic columns						
		- Brackets / mounts						
		-Upholstery and covers						
		- Fastening/clamping elements						
		Pelvi.Loc pelvic belt (only for Benni Light)						
		-Fastening and closure						
		-Textile elements and straps						
5	Back/chest/ head area	Back/pelvic pad						
		- Brackets						
		-Upholstery and covers						
		Fastening/clamping elements						
		Headrests *						
		- Brackets						
		- Upholstery and covers						
		- Fastening/clamping elements						
		The sternum *						
		- Brackets						
		- Upholstery and covers						
		- Fastening/clamping elements						
6	Table area / Armrests *	Mounts						
		Table columns						
		Table mount with depth adjustment						
		Wooden parts						
		Fastening/clamping elements						

(\*) – only for the Heidelberg Benni standing device

Notes for any repairs and further maintenance:


Inspection carried out on

Inspection carried out by

Signature