

Operator

Product

Serial number

Maintenance interval

12 months

Pos.	Assembly		Settings & functions		Damage & deformation		Bolted joints	
			no defects	defects found	no defects	defects found	no defects	defects found
1	Frame	Standard frame, abducted						
		Closed frame, abducted						
		Tip assist						
		Tipping bar						
		Parking brakes/knee lever brakes						
		Brake with extension lever						
		Wheel base extension						
		Castor wheels with wheel forks and castor stems						
		Castor wheel stop						
		Width extension kit for castor wheels						
		Drive wheels with drum brake carrier and hand brake levers						
		Quick-release axles						
		Spoke guards						
		Width extension kit for drive wheels						
		Wheel guards						
		Seat tilting via gas spring						
		One-hand operation of seat tilting						
		Outdoor front wheel accessory incl. wheel chocks						
		Fenders for outdoor front wheel accessory						
		Wheelchair-tiedown and occupant-restraint system						
		Transport fittings						
		scala-/viamobil attachment						
		scalamobil anti-tip						
2	Seat	Seat plate						
		Trapezoid adapter						
		Mounting for footplate						

Pos.	Assembly		Settings & functions		Damage & deformation		Bolted joints	
			no defects	defects found	no defects	defects found	no defects	defects found
		Seating element						
		Back element						
		Connecting hinge						
		Extension of the seating element						
		Universal bracket						
		Headrest, shell-shaped						
3	Back	Back adjustment, mechanical						
		Back adjustment via gas spring						
		Push bar, angle-adjustable						
		Dynamic back						
		Mounting device with a split, fixed bracket						
		Mounting device with a split, adjustable bracket						
		Standard armrests						
		Comfort armrests						
4	Feet	Footplate						
		Light-weight footplate						
		Footplates, parted						
		Foot angle, adjustable						
		Folding mechanism						
		Heel loop						
		Calf support						
		Calf supports, split						
		Flip-up foot rests, standard						
		Flip-up foot rests, special						
		Heel loop for flip-up foot rests						

Notes on repairs and maintenance

Inspection completed on

Inspection completed by

Signature