

# REHATEC®

## Product Overview



We set standards.

Rehatec® GmbH  
In den Kreuzwiesen 35  
69250 Schönau  
Germany

[www.rehatec.com](http://www.rehatec.com)  
[info@rehatec.com](mailto:info@rehatec.com)

©2023 Rehatec® GmbH  
All rights reserved

# Contents

<b>Our website</b> .....	4
--------------------------	---

<b>Our products</b> .....	5
---------------------------	---



## Standing

Heidelberg upright and supine stander Lasse .....	6
---	---

Biometric standing system Da Vinci .....	10
--	----

Heidelberg upright stander Standard .....	12
---	----



Heidelberger Stehgerät Impuls .....	14
-------------------------------------	----

Heidelberg upright stander Jumbo and Jumbo for hip-knee-ankle-foot orthoses .....	16
--	----

Heidelberg upright stander Benni .....	18
--	----

Heidelberg upright stander Speedy .....	20
---	----

## Seating



Pauli activity chair .....	22
----------------------------	----

Nele and Nele Classic activity chair .....	24
--	----

Noah hi-lo base .....	26
-----------------------	----

Nick hi-lo base .....	27
-----------------------	----

## Moving



Galileo mobility base .....	28
-----------------------------	----

Sunny and Sunny Swing mobility base .....	30
---	----

Sunny mobility base suitable for seat shells for supine positions .....	32
---	----

Mika mobility base .....	33
--------------------------	----

## Headrests



Pro-System .....	34
------------------	----




Heads Up .....	35
----------------	----

<b>Colour and fabric overview</b> .....	36
---	----

# Digitally connected

## Website

Our website is full of information such as current order forms, instructions for use and spare parts lists.



# Our products



## Standing

Since we know the importance of standing, developing innovative standing frames is a matter that is very close to our heart. Standing is a major milestone in developing own locomotion skills. By putting users in an upright position, our standing frames improve vital functions, activate muscles, stimulate brain activity, and stabilize bony structures. Users are able to participate in everyday life at eye level. Our upright standers are a perfect fit, meet individual requirements and offer built-in growth options.



Standing frames by Rehatec® accompany users from early childhood to adulthood.



## Seating

Various disease patterns, such as cerebral movement disorder, paralysis, or major developmental delays, require perfect support in a sitting position to improve the user's quality of life, to gain independence and to encourage increased participation.

The expertise of many therapists, specialist dealers and users is always part of our development process since we want to satisfy all therapeutic requirements and ensure a user-friendly, functional application. Our activity chairs are the result of good teamwork.



## Moving

Mobility bases cater to users who are not able to move around on their own and who need a customised seating system. Combined with a customised seat shell that perfectly adapts to the mobility base, users gain a high degree of mobility and benefit from increased participation. This increases the quality of life both at home and outside.

Rehatec® defines all technical requirements e.g., regarding frame, recline or individual adjustments, together with therapists, users and specialist dealers and evaluates the results in close collaboration.



## Headrests

The use of headrests on standing devices, mobility bases, activity chairs or hi-lo bases brings many benefits to patients. When used correctly, they can improve head control, breathing and the ability to focus attention, improve food intake and raise visual and auditory perception. Patients are able to participate more actively in life because it is easier to communicate with others and interact with them.

Rehatec headrests can be individually adapted to patients' needs, from the number of pads and flip away arms to the headband or chin prompt.

---

## Age groups

Infants



to approx.  
3 years old

Teenagers



from approx. 14 to  
18 years old

Children



from approx. 4 to  
13 years old

Adults



from approx.  
18 years old



# Heidelberg upright and supine stander Lasse

## Sizes Mini to 2

The height adjustable upright and supine stander Lasse moves young users from a supine to a standing position and enables them to stand at eye level with their friends. The variable transfer height ensures a safe, stress-free positioning of the user in supine.



Lasse size 1  
Frame colour: pure white  
Fabric: light blue/grey



Lasse size 2  
height adjustable  
Frame colour: anthracite grey  
Fabric: apple green/grey



Degree indicator for height adjustable frame



Leg positioning, angle adjustable



Lateral supports, flip away



Homepage



Product page/  
technical data



Order form  
size Mini



Order form  
size 1 and 2



Colour overview

GMFCS  
Level  
IV - V

3 Sizes

0.65 m -  
1.35 m

Max.  
60 kg

6 different  
covers

22 frame  
colours

**Also available  
with height adjustment**



Infants



Children



Teenagers





# Heidelberg upright and supine stander Lasse

## Sizes 3 to 5

The recline of the Lasse upright and supine stander can be steplessly adjusted and gently brings users of all ages to a standing position. Since the upright and supine stander adapts to users from head to toe, Lasse perfectly meets individual needs.



Lasse size 3  
Frame colour: metallic silver  
Fabric: light blue/grey




Lasse size 4  
Frame colour: metallic silver  
Fabric: apple green/grey



Individual leg abduction



Butterfly harness matches fabric



Knee cups, contoured



Homepage



Product page/  
technical data



Order form  
size 3 and 4



Order form  
size 5



Colour overview

GMFCS  
Level  
IV - V

3 Sizes

1.28 m -  
2.00 m

Max.  
130 kg

10 different  
covers

22 frame  
colours

Also available  
with height adjustment



Teenagers



Adults





# Biometric standing system Da Vinci

The Da Vinci standing system is an upright and supine stander offering a modular design. This design gives even users with severe disabilities the opportunity to stand upright. Each single element is very flexible to adjust and the standing frame perfectly adapts to the user and meets most challenging needs.



Pommel



Adjustable back elements



Headrest, shell-shaped



Homepage



Product page/  
technical data



Order form



Colour overview

GMFCS  
Level  
IV - V

2 Sizes

1.20 m -  
1.80 m

Max.  
90 kg

4 cushion  
colours

22 frame  
colours

Teenagers

Adults





# Heidelberg upright stander Standard

The Standard upright stander is a standing frame for classical standing in an upright position. Combined with the Multi tilting base, it activates even inactive muscle groups. Thanks to its modular design, the stander offers built-in growth options and is easy to adjust.



Standard size 2  
Multi base  
Frame colour: bordeaux violet  
Fabric: alue atoll artificial leather



Standard size 3  
Frame colour: metallic silver  
Fabric: anthracite artificial leather



Tray with bowl and cover



Removable pelvic supports



Multi tilting base



Homepage



Product page/  
technical data



Order form



Colour overview

GMFCS  
Level  
III - IV



3 Sizes



0.65 m -  
1.90 m



Max.  
100 kg



2 different  
covers



22 frame  
colours



Children



Teenagers



Adults





# Heidelberg upright stander Impuls

Thanks to its integrated electric stand aid, the Impuls upright stander moves users slowly, gently and comfortably from a sitting to an upright position. With its broader frame and a low entry, it offers enough space to optimally position a wheelchair and makes the transfer easier and safe. Buttock support belt, supports and other accessories provide stability which is essential for effective standing training.



Frame with low entry



New electric drive with emergency stop button



Back support



Homepage



Product page/  
Technical data



Order form



Colour overview

GMFCS  
Level  
IV



2 Sizes



1.30 m -  
1.95 m



Max.  
130 kg



2 different  
covers



22 frame  
colours



Teenagers



Adults





# Heidelberg upright standers Jumbo & Jumbo for hip-knee-ankle-foot orthosis

The Jumbo upright stander is a standing frame for children's first successful standing training. For children with hip-knee-ankle-foot orthosis, a separate model is available to prevent children from falling. The orthosis is securely attached to the stander.




Jumbo for hip-knee-ankle-foot orthosis  
Frame colour: yellow green



Multi tilting base



Metal loop for hip-knee-ankle-foot orthosis



Tray with elephant design



Homepage



Product page  
Jumbo



Order form  
Jumbo



Product page Jumbo  
für Stehorthesen



Order form Jumbo  
für Stehorthesen

GMFCS  
Level  
III - IV

1 Size

0.65 m -  
1.10 m

Max.  
40 kg

2 different  
covers

22 frame  
colours



Infants



Children






# Heidelberg upright stander Benni

The Benni upright stander provides static and dynamic standing training for children. Thanks to its dynamic centre column, the standing frame teaches them how to maintain their equilibrium and keep their balance. Benni helps children (re)gain their postural stability.



Benni size 2  
Frame colour: heather violet  
Fabric: anthracite



Dynamic centre column



Swivel arm frame



Sternum chest support



Homepage



Product page/  
technical data



Order form



Colour overview

GMFCS  
Level  
III - IV



2 Sizes



0.65 m -  
1.35 m



Max.  
45 kg



2 different  
covers



22 frame  
colours



Infants



Children



Teenagers






# Heidelberg upright stander Speedy

The Speedy upright stander provides children and teenagers with mobility during standing training. Thanks to big drive wheels, users can move around independently. Users apply the mobile stander to train muscles, bones and joints, but also expand their radius of action.



Speedy size 1  
Frame colour: anthracite metallic  
Fabric: blue atoll artificial leather



Chest support with lateral supports



Tilt-in-space adjustment



Knee cups, contoured



Homepage



Product page/  
technical data



Order form



Colour overview

GMFCS  
Level  
III - IV

4 Sizes

0.90 m -  
1.80 m

Max.  
75 kg

2 different  
covers

22 frame  
colours



Infants



Children



Teenagers



# Pauli activity chair

The Pauli activity chair is a child's companion in everyday life and preschool with built-in growth options. Caregivers hardly need any tools to adjust the height, depth or angle of its footplate, seat plate or backrest.

The KiGa version offers a footplate that is particularly low above the ground and it thus adapts perfectly to low table heights.



Pauli size 1 KiGa  
Frame colour: beech  
Fabric: Trevira® grey



Lateral support with protective cover over bracket



Height-/angle-adjustable backrest



Quick height adjustment



Homepage



Product page/  
technical data



Order form



Colour overview

GMFCS  
Level  
I - III



3 Sizes



0.80 m -  
1.50 m



Max.  
50 kg



7 different  
covers



1 frame  
colour



Children



Teenagers



# Nele activity chair


The Nele activity chair that fosters inclusivity and already supports infants in their first attempts at sitting. Thanks to the height adjustment, Nele adapts to various table heights, even the lowest ones. The seat tilting offers a resting position with pressure relief.



Nele Classic size 1  
Frame colour: pure white  
Fabric: blue microfibre



Unique seat depth adjustment



Stepless height adjustment based on the Noah hi-lo base



Powered height adjustment



Homepage



Product page/  
technical data



Order form  
Nele



Order form  
Nele Classic



Colour overview

GMFCS  
Level  
III - IV

4 Sizes

0.65 m -  
1.50 m

Max.  
60 kg

7 different  
covers

22 frame  
colours

Size Mini now  
available for toddlers



Infants



Children



Teenagers



Order form  
Nele Mini



Order form  
Nele Classic Mini



# Noah hi-lo base

Infants



Children



Teenagers



The Noah hi-lo base is a tilting frame for all common seat shells. Noah offers gas spring, hydraulic or powered height adjustment for an optimal seat height at all times. With outdoor tyres, it is also suitable for outdoor activities.



Noah size 2  
Frame colour: ultramarine blue

GMFCS  
Level  
IV - V



2 Sizes



Max.  
75 kg



22 frame  
colours



Outdoor frame



Height adjustment to  
lower the seat down to the  
ground



Trapezoid adapter



Homepage



Product page/  
technical data



Order form



Colour overview



# Nick hi-lo base

The Nick hi-lo base can be combined with all common seat shells. Thanks to its manual height adjustment, Nick can be adapted to many table heights. The seat tilt is steplessly adjustable by means of gas spring.



Infants



Children



Teenagers



Nick size 3  
Frame colour: ultramarine blue


GMFCS  
Level  
IV - V

3 Sizes



Max.  
90 kg

22 frame  
colours



Adjustable back attachment



Steplessly adjustable seat tilt



Padded footplate



Homepage



Product page/  
technical data



Order form



Colour overview



# Galileo mobility base

The versatile Galileo mobility base is suitable for all common seat shells. It offers adjustable seat tilt, dynamic back adjustment and a displaceable seat plate to ensure an optimal sitting position at all times. The extensive range of accessories meets almost any requirements.



Flip-up leg positioning



Vent tray



Wheelchair-tiedown and occupant-restraint system



Homepage



Product page/  
technical data



Order form  
Galileo



Colour overview

GMFCS  
Level  
III - V

8 Sizes

Max.  
120 kg

22 frame  
colours



Children



Teenagers



Adults





## Sunny & Sunny Swing mobility bases

The Sunny or Sunny Swing mobility bases are suitable for all common seat shells. Due to variable options to install the wheels, they are even suitable for extremely narrow or wide seat shells. The back is angle adjustable and can be folded forwards or backwards.




Sunny size 2  
Abduction frame  
Frame colour: pure orange



Parking/knee lever brake



Anti-tip



Front castor with fork

GMFCS  
Level  
IV - V

4 Sizes

Max.  
120 kg

22 frame  
colours

Children

Teenagers

Adults

Sunny Swing size 3  
Frame colour: ruby red



Homepage



Product page/  
technical data



Order form



Colour overview




# Sunny mobility base

Suitable for seat shells for supine positions


Thanks to its sturdy frame and extendable wheelbase, the Sunny mobility base is ideal for combination with individual seat shells for supine positions. The angle adjustable back and footrests give users the opportunity to take an individual pressure-free lying position



Sunny for supine positions size 3  
Frame colour: graphite black matt



Backrest, angle adjustable



Footrest, angle adjustable



Homepage



Product page/  
technical data



Order form  
Sunny for supine positions



# Mika mobility base

The sturdy Mika mobility base can be combined with all common seat shells. The wheelbase is adjustable in six positions to offer different driving characteristics - from sporty to extremely stable. The back is angle adjustable and can be folded down.



Mika size 2  
Frame colour: pure white



Children



Teenagers



Adults

GMFCS  
Level  
IV - V



4 sizes



Max.  
120 kg



22 frame  
colours



Seat tilt



Tip assist



Parking/knee lever brake




Homepage



Product page/  
technical data



Order form



Colour overview




## Pro-System headrests




Whitmyer's versatile Pro-System range of headrests consist of 1 to 4 modules. The modules are steplessly adjustable and rotatable for optimal adaptation. Their curved shape supports the back of the head and does not restrict acoustic perception.



Headrest, adjustable



Headrest with integrated lateral wings



Headrests with adjustable lateral wings



Homepage



Product page/  
technical data



Order form  
1-/2-modules



Order form  
3-/4-modules



# Heads Up headrests


The Heads Up headrest system features flip away, adjustable arms. This means that the ears are not touched when positioning or moving users. The headrest is available in different sizes, collar sizes and arm lengths.



Chin prompt



Occipital support



Flip away arms



Homepage



Product page/  
technical data




















Order form

# Fabrics, fabric colours and frame colours

We offer different fabrics for different needs.

All frames are powder coated. The colours are based on the standardised RAL system.

	Microfibre <b>Lilac</b>		Microfibre <b>Blue</b>		Microfibre <b>Green</b>		Trevira® CS <b>Yellow</b>
	Trevira® CS <b>Pink</b>		Trevira® CS <b>Red</b>		Trevira® CS <b>Apple green</b>		Trevira® CS <b>Light blue</b>
	Trevira® CS <b>Grey</b>		Trevira® CS <b>Black</b>		Artificial leather <b>Sand</b>		Artificial leather <b>Terracotta</b>
	Artificial leather <b>Dark red</b>		Artificial leather <b>Dark blue</b>		Artificial leather <b>Light grey</b>		Artificial leather <b>Anthracite</b>
	Artificial leather <b>Black</b>	You can find the available fabrics for your desired device on our homepage or in the respective order form. Other colour combinations are available on request subject to an extra charge.					

	<b>Pure white</b> RAL 9010		<b>Daffodil yellow</b> RAL 1007		<b>Luminous yellow</b> RAL 1026		<b>Pure orange</b> RAL 2004
	<b>Luminous orange</b> RAL 2005		<b>Ruby red</b> RAL 3003		<b>Pearl ruby red</b> RAL 3032		<b>Heather violet</b> RAL 4003
	<b>Bordeaux violet</b> RAL 4004		<b>Ultramarine blue</b> RAL 5002		<b>Metallic traffic blue*</b>		<b>Pearl night blue</b> RAL 5026
	<b>Water blue</b> RAL 5021		<b>Light green</b> RAL 6027		<b>Yellow green</b> RAL 6018		<b>Pearl green</b> RAL 6035
	<b>Light green</b> RAL 6038		<b>Anthracite grey matt</b> RAL 7016m		<b>Metallic silver*</b>		<b>Metallic anthracite*</b>
	<b>Graphite black matt</b> RAL 9011m		<b>Glossy jet black</b> RAL 9005g				

\*no RAL colours

