

# Heidelberg Upright Stander Standard **REHATEC®**

suitable from about 100 cm to about 190 cm with dynamic functions

## Cost estimate

Sender / billing address



## Order

Address for dispatch

Purchaser

Date

Commission



### Base Model Size 2/3 (102/101)

- Dynamic functions (movable footplate, movable/dynamic knee cups, movable swivel arm frame)
- Castors 75 mm with lock
- Sandals, screw-fastening, half-height
- Therapy table with a semicircular mould (standard)
- Therapy table adjustable for depth

#### Intended age group:

For size 2, approx. 3 to approx. 14 years of age\*

for size 3, approx. 14 years\*

**Max. user weight:** 100 kg

**User height:** for Size 2 about 100 cm to 140 cm  
for Size 3 about 135 cm to 190 cm

Figure does not comply with the base model and its accessories | Size tolerance of  $\pm 3\%$

\*Appropriate patient measurements are crucial for patient care.

Additional technical details available at [www.rehatec.com](http://www.rehatec.com) -> Downloads. Or send us your request by email, fax or letter.


Valid from 01.01.2026 until recalled. © Unauthorized duplication without prior permission of Rehatec GmbH is strictly prohibited.

Technical changes reserved. Our General Terms and Conditions apply, see [www.rehatec.com](http://www.rehatec.com).

By ordering our products, you confirm that you have read and accepted our General Terms and Conditions

Rehatec® GmbH, In den Kreuzwiesen 35, 69250 Schönau, T: +49 6228 9136-0, F: +49 6228 9136-99, [info@rehatec.com](mailto:info@rehatec.com), [www.rehatec.com](http://www.rehatec.com)

Please enter the measurements of your patient.



**\* CAUTION:**  
For technical reasons, the distance between the tray height (a) and the height of the pelvic middle (c) amounts to at least 15 cm!

Size tolerances ±3%

Base model/user measurements					
Size	Tray height (a)	Knee height (b)	Height pelvic middle (c)	Knee width (e)	Pelvic width (f)
2	70 - 90 cm	20 - 45 cm	50 - 70 cm	6 - 14 cm	26 - 43 cm
3	85 - 119 cm	27 - 56 cm	70 - 104 cm	6 - 14 cm	29 - 53 cm
User	cm	cm	cm	cm	cm

Size	Chest width (g)	User weight
2	20 - 32 cm	max. 100 kg
3	30 - 42 cm	max. 100 kg
User	cm	kg

Equipment measurements			
Size	Total length	Inner frame width	Total width
2	100 cm	49 cm	61 cm
3	115 cm	60 cm	75 cm <sup>1</sup>

<sup>1</sup> Selecting the broader frame, the total width amounts to 84 cm!

### Base model

**Heidelberg Standing Frame Standard size 2/3**

Size	Code
2 (102)	1017707-1
3 (101)	1017350-1

<b>Size 2</b>	<b>6.724,-€</b>
<b>Size 3</b>	<b>6.724,-€</b>

### Dynamic functions

**Dynamic functions incl. movable footplate, movable/dynamic knee cups, movable swivel arm frame**

<b>Size 2</b> (Code 1026699)	<b>Standard</b>
<b>Size 3</b> (Code 1026102)	<b>Standard</b>

### Frame colours

Pure white (RAL 9010)	Pearl night blue, matt (RAL 5026)
Daffodil yellow (RAL 1007)	Water blue (RAL 5021)
Luminous yellow (RAL 1026)	Light green (RAL 6027)
Pure orange (RAL 2004)	Yellow green (RAL 6018)
Luminous orange (RAL 2005)	Pearl green, matt (RAL 6035)
Ruby red (RAL 3003)	Luminous green (RAL 6038)
Pearl ruby red, matt (RAL 3032)	Anthracite grey, matt (RAL 7016m)
Heather violet (RAL 4003)	Metallic silver* <b>*no RAL colours</b>
Claret violet (RAL 4004)	Metallic anthracite*
Ultramarine blue (RAL 5002)	Graphite black, matt (RAL 9011m)
Traffic blue (RAL 5017)	Glossy jet black (RAL 9005g)

**This is only an overview of colours. Find all details at [www.rehatec.com](http://www.rehatec.com).**

### Plain cushion colours

- Light grey artificial leather (N20057)
- Anthracite artificial leather (N20058)

**Additional faux leather colours available on request for an additional charge of €240 (special construction).**

Frame		Wall and furniture bumpers for 75 mm castors	Frame, broad
<b>Lockable castors, 75 mm</b> Code 1005593-1 	<b>Wall and furniture bumpers for 75 mm castors</b> Code 1018607 	<b>Frame, broad</b> Code 1027159 incl. footplate, broad Internal dimension (a) with 75 mm castors = 68 cm, total width = 84 cm 	<p><b>Using the broader base frame, the total width amounts to 82 cm</b></p> <p><b>Only for Upright Stander Standard Size 3!</b></p>
<b>Standard</b>	<b>Set</b>	<b>125,-€</b>	<b>Extra charge</b>
			<b>614,-€</b>




Additional technical details available at [www.rehatec.com](http://www.rehatec.com) -> Downloads. Or send us your request by email, fax or letter.



Valid from 01.01.2026 until recalled. © Unauthorized duplication without prior permission of Rehatec GmbH is strictly prohibited.

Technical changes reserved. Our General Terms and Conditions apply, see [www.rehatec.com](http://www.rehatec.com).



By ordering our products, you confirm that you have read and accepted our General Terms and Conditions

Rehatec® GmbH, In den Kreuzwiesen 35, 69250 Schönau, T: +49 6228 9136-0, F: +49 6228 9136-99, info@rehatec.com, www.rehatec.com



Therapy table		Depth adjustable therapy table	Table skirting											
<b>Standard therapy table</b>  <table border="1"> <tr><th>For Size</th><th>Code</th></tr> <tr><td>2</td><td>1004331</td></tr> <tr><td>3</td><td>1017349</td></tr> </table>	For Size	Code	2	1004331	3	1017349	<b>Depth adjustable therapy table</b> Code 1018797-1 	<b>Table skirting</b>  <table border="1"> <tr><th>For Size</th><th>Code</th></tr> <tr><td>2</td><td>1018800</td></tr> <tr><td>3</td><td>1018799</td></tr> </table>	For Size	Code	2	1018800	3	1018799
For Size	Code													
2	1004331													
3	1017349													
For Size	Code													
2	1018800													
3	1018799													
<b>Standard</b>	<b>Standard</b>	<b>Satz</b>	<b>173,-€</b>											

Tray height adjustment by gas spring					Handles (vertical)						
<b>NEW</b>					<b>NEW</b>						
					Code 1027066 						
For Size	Code	Tray height (a) <sup>1</sup>	Knee height (b)	Height pelvic middle (c)	<table border="1"> <tr><th>Code</th><th>Montage</th></tr> <tr><td>1027066-1</td><td>enclosed separately for self-assembly</td></tr> <tr><td>1027066-2</td><td>Install grab rails</td></tr> </table>	Code	Montage	1027066-1	enclosed separately for self-assembly	1027066-2	Install grab rails
Code	Montage										
1027066-1	enclosed separately for self-assembly										
1027066-2	Install grab rails										
2	1026107-A	ca. 81 – 90 cm	20 – 40 cm	63 – 73 cm							
3	1025994-A	ca. 91 – 115 cm	28 – 60 cm	72 – 100 cm							
	patient measurements <sup>2</sup>	cm	cm	cm							
<b>Set (customized)</b>					<b>1.009,-€</b>						
<b>Pair</b>					<b>195,-€</b>						

<sup>1</sup>Only possible with swivel arm basin frame. <sup>2</sup>Automatically transferred from page 2.

Sandals		Foot strap																					
<b>Sandals, screw-fastening, half-height</b>  <table border="1"> <tr><th>Size</th><th>Code</th><th>a</th><th>b</th><th>c</th></tr> <tr><td>1</td><td>1022897</td><td>11 cm</td><td>14 cm</td><td>3 cm</td></tr> <tr><td>2</td><td>1022898</td><td>18,5 cm</td><td>19 cm</td><td>4 cm</td></tr> </table>		Size	Code	a	b	c	1	1022897	11 cm	14 cm	3 cm	2	1022898	18,5 cm	19 cm	4 cm	 <table border="1"> <tr><th>Size</th><th>Code</th></tr> <tr><td>Short version</td><td>1014394</td></tr> <tr><td>Long version</td><td>1014395</td></tr> </table>	Size	Code	Short version	1014394	Long version	1014395
Size	Code	a	b	c																			
1	1022897	11 cm	14 cm	3 cm																			
2	1022898	18,5 cm	19 cm	4 cm																			
Size	Code																						
Short version	1014394																						
Long version	1014395																						
<b>Standard</b>		<b>Pair</b>																					
		<b>61,-€</b>																					

**Recommended for screw-on foot shells**

Knee cups		Neoprene cover, black, for movable knee cups
<b>Knee cups made to measure</b> Code 10045-2 	<b>Neoprene cover, black, for movable knee cups</b> Code 1026799 	<b>NEW</b>
Size e = _____ cm (from 15 cm)	<b>Size complies with user measurements (6 – 14 cm)</b>	
<b>Pair / extra charge (customized)</b>	<b>216,-€</b>	<b>Pair / extra charge</b>
		<b>175,-€</b>

Additional technical details available at [www.rehatec.com](http://www.rehatec.com) -> Downloads. Or send us your request by email, fax or letter.

Valid from 01.01.2026 until recalled. © Unauthorized duplication without prior permission of Rehatec GmbH is strictly prohibited.

Technical changes reserved. Our General Terms and Conditions apply, see [www.rehatec.com](http://www.rehatec.com).

By ordering our products, you confirm that you have read and accepted our General Terms and Conditions

Rehatec® GmbH, In den Kreuzwiesen 35, 69250 Schönau, T: +49 6228 9136-0, F: +49 6228 9136-99, info@rehatec.com, www.rehatec.com

